Servotome

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ENGLISH

C O N T E N T S

I-	BASIC FEATURES	6
II-	INSTALLATION	8
III-	RECOMMENDATIONS FOR USE	8
IV-	USE	11
V-	MAINTENANCE	12
VI-	TROUBLESHOOTING	13
VII-	TECHNICAL CHARACTERISTICS	15

I - BASIC FEATURES

The SERVOTOME® is a high-frequency current generator used for the incision and coagulation of soft gingival tissue thanks to a wide range of electrodes.

1. 1 GENERAL CHARACTERISTICS

- Incision by high-frequency current
- Coagulation by pulsating current
- Status displayed by indicator lamps plus audible signal
- Power control button
- Coagulation control button
- Control by footswitch.

1. 2 ITEMS DELIVERED

- One SERVOTOME® device
- One footswitch and lead
- One removable mains lead with earthed connection
- One electrode-holder
- One electrode-holder cord
- One capacitive coupler (back plate)
- One box containing 10 electrodes
- One user's manual
- Setting instructions.

1. 3 GENERAL VIEW

Refer to the illustrations at the end of this manual.

1	Rear Panel	6	Electrode-holder
2	Indicator Lamps	7	Electrode-holder rest
3	Power Control Button	8	Electrode-holder cord
4	Coagulation Control Button	9	Footswitch
5	Capacitive coupler		

1. 4 FRONT PANEL

This comprises two indicator lamps and two control buttons.

- "Active" Indicator Lamp (10)

This yellow indicator remains lit as long as the footswitch is pressed and indicates the presence of the high-frequency current. Note that a buzzer sounds while the footswitch is being pressed.

- ON Indicator Lamp (11)

This is green; it lights up when the ON/OFF switch (16) of the device is ON.

- Button  (12)

Minimum to maximum incision power setting. At maximum power (setting 10), the power delivered is 50W; this does depend, however, on the operating conditions and the histological variables.

- Button (13)

Coagulation Control:

- Value 1: Minimum coagulation
- Value 10: Maximum coagulation.

1. 5 REAR PANEL

- Capacitive coupler Connector (14)

Connection to the capacitive coupler: please refer to "Recommendations for use" for the position of the capacitive coupler.

- Footswitch Connector (15)

Connection to the single footswitch.

- ON/OFF Switch (16)

Position I: ON

- Fuse Housing (17)

Contains the fuses designed to protect the unit.

- Mains Plug (18)

Connection to the mains by a standard plug.

1. 6 SINGLE FOOTSWITCH

Pressing the footswitch (9) activates the high-frequency output.

1. 7 CAPACITIVE COUPLER

The capacitive (5) capacitive coupler connects to the device by a lead and returns the high-frequency current. It should be placed on the clothed patient's back (scapula).

II - INSTALLATION

2. 1 RECOMMENDATIONS

For your safety, connection to the electrical circuit should be carried out by an approved dental installation engineer.

Important: before using the device for the first time, clean, decontaminate and sterilize the accessories (in particular the electrode-holder and the electrodes). See Section 5 - Maintenance.

2. 2 SETTING UP THE EQUIPMENT

- Place the device where it is to be used, on a level surface or on one that has a slope of less than 10°.
- Plug in the electrode-carrier cord (8).
- Connect the footswitch (9) to the footswitch connector (15).
- Connect the capacitive coupler (5) to its socket (14).
- Connect the mains plug into the mains socket (18). It is essential for the mains socket to be earthed.

III - RECOMMENDATIONS FOR USE

3. 1 GENERAL RECOMMENDATIONS

Before the SERVOTOME® is used for the first time, it is essential to carry out tests on anatomical parts (pieces of meat - ideally, a piece of beef heart, chicken breast, etc.) to determine how they react to incision and to help adopt the right clinical procedure (electrode movement speed). Repeat these exercises as many times as necessary.

3. 2 THE DEVICE

Make sure the device is properly plugged in and insulated. Use of the mains socket with its earthed connection is mandatory.

The user should check the accessories and leads regularly so as to locate any insulation fault and replace them if and when necessary.

Do not use the device if it appears to be defective in any way.

In the event of a malfunction, you are advised to contact your SERVOTOME® supplier rather than an official repairs technician.

3. 3 ELECTRODES AND ELECTRODE-HOLDERS

It is essential to replace used or damaged SERVOTOME® accessories with the appropriate new SATELEC accessories.

To avoid damaging or burning out the SERVOTOME®, it should not be used with any accessories other than those supplied by SATELEC.

Prior to each use, the electrodes and electrode-holder must be cleaned, decontaminated and sterilized.

The electrode-holder cord must be placed so as to avoid any contact with the patient or other conductors.

Active electrodes that are not being used temporarily must be kept well away from the patient.

The patient should never come into contact with metal parts connected to the earth (operating table, support, etc.).

Contact between the practitioner's skin and that of the patient should be avoided.

Simultaneous use on a patient of a monitoring device equipped with needle electrodes and an electro-surgical device is not recommended; if the need arises, the electrodes in the monitoring device must be placed as far as possible from the electro-surge electrode.

Monitoring systems that include high-frequency current limiters are recommended.

During surgical procedures in which the high-frequency current could flow through parts of the body with a relatively small straight surface section, the use of bipolar methods may be desirable in order to avoid unwanted coagulation.

3. 4 CAPACITIVE COUPLER

Make sure that the entire surface of the capacitive coupler is in contact with the patient's clothed body. The capacitive coupler should be placed between the dental chair and the clothed patient.

The entire surface of the capacitive coupler should be connected as near as possible to the operating area, so that the current's shortest path between the active electrode and the capacitive coupler is not via the heart.

3. 5 OUTPUT POWER

The output power must be the lowest possible for the purpose required.

Abnormally low power can be caused by a bad capacitive coupler connection. Before increasing the power, the operator should check the connections of the different electrodes.

The power level delivered by the SERVOTOME® may vary according to numerous parameters (operating speed, histological variables, thickness and material of clothing, etc.).

A malfunction may cause an abnormal increase in the output power.

3. 6 ENVIRONMENT

Due to the high-frequency rays, there may be some disturbance caused to other electrical equipment when the SERVOTOME® is in use.

The SERVOTOME® is not designed to function in the presence of an anesthetic gas.

The operator must be made aware of the risks involved in the use of high-frequency equipment when flammable gases are used (explosion hazard).

For patients and/or users carrying a pacemaker or any other active implant, there is a potential risk due to possible interference affecting the pacemaker or damage to it. If there is any doubt, seek the advice or approval of suitably qualified persons.

Certain body cavities and/or certain oxygen-charged materials such as cotton, wool and gauze can be ignited by sparks produced during normal use of the SERVOTOME®.

3. 7 CONTRAINDICATIONS

The SERVOTOME® should not be used:

- in the presence of patients or practitioners carrying active implants such as a pacemaker without the advice or approval of suitably qualified persons.
- in the presence of a flammable anesthetic gas.
- with any accessories other than those supplied by SATELEC.

Do not check to see if you have a high-frequency current by producing electric arcs on metal components.

The SERVOTOME® is not designed to withstand defibrillator discharges.

The user should not attempt to repair or modify the unit without prior SATELEC authorization.

3. 8 ADDITIONAL INFORMATION

This information is intended to acquaint the operator with the limitations of electro-surgery.

- Unruly, emotional or excessively nervous patients
- Incomplete anesthesia
- Delicate surgery (mucoperiosteal surgery, grafts, etc.)
- Very fragile tissue
- Ignorance of the theory of electro-surgery
- Lack of practice on anatomical parts
- Ignorance of the patient and of his/her general condition

3. 9 LIABILITY

The manufacturer may not be held liable in any of the following cases:

- Failure to follow the manufacturer's recommendations at the time of installation (mains voltage)
- Work including repairs carried out by persons not approved by the manufacturer
- Use of electrical equipment that does not comply with the regulations in force
- Any use other than the uses specified in this manual
- Failure to comply with the requirements specified in this manual

The manufacturer reserves the right to modify the unit and/or amend the user manual without prior notice.

IV - USE

Use of the SERVOTOME® is restricted to capable qualified practitioners (dental surgeons) and then only in a dental surgery.

Prior to initial use:

- Clean, decontaminate and sterilize the accessories (in particular the electrode-holder and the electrodes). See section entitled Maintenance: Sterilization.
- Perform tests on pieces of anatomy (pieces of meat - ideally, a piece of cow heart, chicken breast, etc.).

4. 1 MOUNTING THE ELECTRODE

- 1- Unscrew the electrode-holder cap (21) a few turns.
- 2- Insert the appropriate electrode (19) for the operating procedure.

For the choice of electrode, refer to the document entitled Servotome® Clinical Instructions.

Important: Do not use an electrode if its plastic sheath looks damaged (splits, holes, etc.) or is missing. Replace the electrode. It is essential to push the electrode well in so that no metal part is visible between the electrode-holder cap and the electrode's plastic sheath (20). Any visible part would cause the current to flow and result in

a painful incision in the wrong part of the patient's mouth.

Replace the electrode-holder if it no longer holds the electrodes tightly.

- 3- Screw the cap back on again (21).

4. 2 START-UP

- 1- Place the capacitive coupler (5) between the chair and the patient's back (scapula). Follow the indication "Patient Side" printed on the capacitive coupler (22). Check to make sure the capacitive coupler is in position on the patient's scapula.

Important: the capacitive coupler does not need to be in direct contact with the patient's skin. It is essential that you refer to the section entitled "Recommendations" before proceeding.

- 2- Turn the mains switch (16) to the ON position (I). The green indicator lamp (11) opposite lights up. The device is on and ready for use.
- 3- Adjust the incision power (▲ button, 12) and coagulation (■ button, 13) settings before operating, to avoid a burn or other undesirable effect.

Important: It is essential that you refer to the section entitled "Recommendations" before proceeding.

ELECTRODE			SECTION	SECTION + MIN. TO MAX. COAGULATION						COAG.	FULG.
			■ ■ ■ =1	■ ■ ■ =2	■ ■ ■ =4	■ ■ ■ =6	■ ■ ■ =8	■ ■ ■ =10	■ ■ ■ =1	■ ■ ■ =1	
I22S	I22CA		▲ =3	▲ =3/4	▲ =3/4	▲ =4/5	▲ =5/6	▲ =5/6			
TR22T	TR22R	TR22L	▲ =4	▲ =4/5	▲ =4/5	▲ =5	▲ =5/6	▲ =6			
I40S	I40CA		▲ =3/4	▲ =3/4	▲ =4	▲ =4	▲ =5/6	▲ =5/6			
FC10N										▲ =5/6	
FC25B	FC32B								▲ =5	▲ =6/7	
▲ = Power (values provided for information only)							■ ■ ■ = Adjustment of coagulation (values provided for information only)				

V - MAINTENANCE

- 4- Press the footswitch (9). Incision is activated immediately. The orange indicator lamp (10) lights up and the buzzer sounds. The indicator lamp goes out when the footswitch is released.

Important: The device should not be used for a period exceeding 10 seconds, following which it should remain inactive for 30 seconds before the operation can be continued.

4. 3 COMPLETE SHUT-OFF

Turn the ON/OFF switch (16) to 0.

5. 1 CLEANING THE CASING

It is not advised to use cleaning products or disinfectants that contain flammable agents. If they do however, make sure the product evaporates and that no traces of the combustible agent remain on the device or its accessories (electrode-holder, electrodes, leads, etc.) when starting it up.

5. 2 ELECTRODES AND ELECTRODE-HOLDER

a) Cleaning the Electrodes

Use alcohol or disinfectant towelettes. Remove any fragments of tissue that may have stuck to the electrodes.

b) Cleaning the Electrode-holder and the Cord

These can be cleaned and disinfected using alcohol or disinfectant towelettes. Make sure the product has evaporated completely. The hollows in the screw thread can be cleaned with cotton wool soaked in alcohol.

Do not attempt to clean the electrode-holder in an ultrasonic tank.

c) Decontaminating the Electrodes

Soak the electrodes in a glutaraldehyde or quaternary ammonium alkaline solution for approximately 20 minutes depending on the concentration (see manufacturer's instructions).

5. 2 STERILIZATION

The cord connecting the electrode-holder to the SERVOTOME® cannot be sterilized.

It is essential to sterilize the electrodes and the electrode-holder before each use.

VI - TROUBLESHOOTING

Requirements and sterilization cycle according to the equipment:

- Equipment: autoclave
- Wet heat sterilization: 134°C at 2 bar (200 KPa) for 18 minutes

Important: metal parts of different types should not touch one another during sterilization. Any such contact would produce electrolytic couples that would cause local damage. To avoid this phenomenon, place the objects in a cloth or a micro-porous sterilizing bag.

Please refer to the table on the following page

5. 3 FUSE REPLACEMENT

Proceed as follows:

- 1- Turn the ON/OFF button to "0" (23) and unplug the mains lead from its socket (24).
- 2- Insert the tip of a screwdriver into the notch (25) to open the protective cap.
- 3- Replace the fuses (26) with fuses of the same rating and type.
- 4- Push the fuse compartment back into its housing.

The SERVOTOME® designed by SATELEC is a medical device of class IIb (European Directive no. 93/42/EEC concerning medical devices; Transposition to the French national legislation by Decree no. 95/292 of 16/03/1995). By application of the French Decree of 5 December 2001 and the corresponding Order of 3 March 2003 concerning the obligation to ensure proper maintenance and quality control of medical devices, the operator -who is responsible for due implementation of the applicable maintenance operations, is invited to refer to, and to apply the maintenance operations as introduced in the following standards:

CEI 1289-1: High-frequency surgical equipment - Part 1: Operation

CEI 1289-2: High-frequency surgical equipment - Part 2: Maintenance.

FAULT	POSSIBLE CAUSE	CORRECTIVE ACTION
Nothing working; the green indicator lamp doesn't light up	The ON/OFF switch is on 0	Turn the ON/OFF switch to I
	Mains lead not properly connected	Check the mains plug
	No electric current	Check electric current supply
	Fuse(s) not working	Change the external fuses. <i>Caution:</i> the unit includes a third, internal, fuse that is inaccessible to the user. Contact the dealer
The green indicator lamp lights up but the orange indicator lamp does not, and the buzzer does not work	Footswitch lead not properly connected	Check the footswitch: plug in firmly
	Faulty footswitch	Replace the footswitch or contact the dealer
	Thermal protection of device	Wait for the device to cool down
The indicator lamps come on and the buzzer is working, but there is no high-frequency current	Electrode-holder not properly connected	Check the electrode-holder cord: - Push unit-side connector in properly - Push electrode-holder-side connector in properly
	Capacitive coupler not properly connected	Check the capacitive coupler: the connector must be all the way in
	Other possibilities	Contact the dealer
The electrode incises with difficulty or is not incising at all	Capacitive coupler in the wrong position.	Follow the indication <i>Patient Side</i> printed on the capacitive coupler. Make sure the entire surface of the capacitive coupler is in contact with the lightly-clothed patient, for example, on the scapula.
	Soiled electrode	Set the unit to 0. Clean the electrode.
	Electrode moving too fast	Reduce operating speed
	Unsuitable electrode	Select the appropriate electrode for the operating procedure
	Electrode worn	Replace the electrode
The electrode sticks to the tissue	Power set too low	Increase power up to the incision threshold. There is no need to go beyond this threshold
	Inappropriate electrode	Select the appropriate electrode for the operating procedure
The electrode incises but there are sparks	Incision power set too high	Reduce incision power down to the threshold. There is no need to go beyond this threshold

VII - TECHNICAL CHARACTERISTICS

7. 1 INPUT POWER

230 V AC $\pm 10\%$ or 115 V AC $\pm 10\%$
Frequency: 50 to 60Hz
Power Consumption: 230 VA

7. 2 OUTPUT POWER

See output power curves below.

7. 3 WORKING CHARACTERISTICS

Monopolar mode bi-terminal technology.
Output power: 50W max. at 375Ω impedance.
Useful output impedance range: 50Ω to $2K\Omega$.
Output frequency: $1.8\text{ MHz} \pm 0.7\text{ MHz}$.
Appliance Class 1 Type BF.
Incision and coagulation settings: 1 to 10.
Maximum output voltage: 725 V peak-to-peak @
 $P = 10, W = 1$.

7. 4 PROTECTION

Two 5mm x 20mm time-delay fuses:
2 x 1.25 AT for the 230 V AC model or
2 x 2 AT for the 115 V AC model.
One 5mm x 20mm internal time-delay fuse: 500 mAT
ref. F1 on the printed circuit (not user-accessible).

Thermal protection in the event of intensive use
and/or non-compliance with the working cycle (10
seconds ON / 30 seconds OFF).

7. 5 USER INTERFACE

ON/OFF switch on rear panel.
Output power control button, settings 1 to 10.
Coagulation control button, settings 1 to 10.
Tight footswitch control switch.
Green ON indicator lamp.
Yellow indicator lamp and buzzer (65dB min. at 1
meter) of HF signal.

7. 6 CONNECTORS

6.35mm capacitive coupler jack connector.
4mm electrode-holder cord safety connector.

7. 7 DIMENSIONS AND WORKING CHARACTERISTICS

Overall dimensions:

- Width: 250mm
- Height: 110mm
- Depth: 240mm.

Weight: 1.2kg.

Fixing: table unit on a level surface, slope not to
exceed 10° .

Working temperature : 10 to 40°C

Working relative humidity: 30 to 75%

Storage/transport temperature: - 40 to $+70^\circ\text{C}$

Storage relative humidity: 10 to 95%.

7. 8 SYMBOLS AND MARKING

Symbols	Meaning
	Alternating current
	BF type appliance
	Caution! Refer to accompanying documents
	Capacitive coupler connector
	Footswitch connector
	Not to be disposed of as household waste
	Earth-referred capacitive coupler for HF currents
CE 0459	CE marking in accordance with Directive 93/42/EEC
230V : 2x1.25 AT 115V : 2x2 AT	Primary fuse rating depending on mains supply
	Power setting control
	Coagulation setting control
IP21	Degree of protection afforded by the casing

7. 9 STANDARDS

Manufactured in accordance with the regulations and standards in force.

European Directive CE 93/42/EEC

IEC 601-1.

IEC 601-2-2.

7. 10 DISPOSAL AND RECYCLING

The SERVOTOME® and its accessories can only be disposed of by an authorized collector of electronic devices or by returning them to the SATELEC After Sales Department.

Under no conditions should the Servotome and its accessories be disposed of by the user.

NB: SATELEC will provide the technical staff of authorized SATELEC dealers upon request with all the necessary information to carry out repairs to parts of the device that SATELEC deems repairable.



The use of the WEEE symbol  indicates that this product must not be treated as household waste, but must be treated separately, in conformity to the Directive 2002/96/CE. By ensuring this product is disposed of correctly, you will help to protect the environment. For more detailed information about the recycling of this product, please contact your local authority, your house waste disposal service provider or where you purchased the product. At the end of their life, tips should be discarded in safety containers designed for this purpose

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