Extraction

User's manual • Notice d'utilisation

PIEZOTOME™ - PIEZOTOME™ 2

IMPLANT CENTER - IMPLANT CENTER 2

ENGLISH

CONTENTS

IMPORTANT : Surgical tips designed for PIEZOTOME and IMPLANT CENTER cannot be used with PIEZOTOME 2 and IMPLANT CENTER 2 and vice versa.

I	Introduction	5
II	Description and clinical applications	6-7
III	Method of operation	7-8
IV	Sterilization	8-10
V	Replacing tips	11
	Subsidiaries	20-22

I - INTRODUCTION

Specially designed to address periodontal and pre-implant surgery, the EXTRACTION Kit reduces post-extraction trauma.

This kit is intended for clinical acts such as extraction of wisdom teeth, fractured roots, impacted teeth and root section. The six ultrasonic tips are each adapted to the anatomical context (reduced visibility and space) of anterior and posterior areas.

The tip fits between the root and the alveolar bone, detaching the ligament. The tooth is then separated and removed "gently". Micro-ultrasonic oscillations cause detachment of the ligament without damaging the surrounding soft tissue.

The sterile spray helps to cool the tips (thus avoiding tissue damage from exposure to high temperatures) and provides excellent visibility of the operating field.

⚠ The PIEZOTOME tips and PIEZOTOME handpiece are designed for use with the PIEZOTOME generator or IMPLANT CENTER combined generator and bone surgery micro-motor, manufactured by Satelec.

The PIEZOTOME 2 tips and PIEZOTOME 2 handpiece are designed for use with the PIEZOTOME 2 generator or IMPLANT CENTER 2 combined generator and bone surgery micro-motor, manufactured by Satelec.

II - DESCRIPTION AND CLINICAL APPLICATIONS

2.1 LC1 TIP

This tip allows cautious and in-depth penetration, along the desmodontium, between the root and the alveolar bone.

2.2 LC1-90° TIP

This ultrasonic periotome facilitates the enlargement of the desmodontal space in interproximal, lingual and distal areas of molars.

2.3 LC2 TIP

The remarkably slim LC2 provides access to reduced spaces between the root and the alveolar bone without damaging the cortical bone.

2.4 LC2L TIP

45° left oriented (view from below), the LC2L facilitates access to posterior areas. The tip's slim line respects the morphology of the table bone and makes it easier to work in areas of reduced visibility.

2.5 LC2R TIP

45° right oriented (view from below), the LC2R facilitates access to posterior areas. The tip's slim line respects the

morphology of the table bone and makes it easier to work in areas of reduced visibility.

2.6 NINJA™ TIP

This double serrated saw with laser markings every 3mm, to evaluate the depth or thickness of cut, is particularly effective for hemisections and root

sections. It is also indicated for sectioning impacted molars during extraction.

III - METHOD OF OPERATION

3. 1 FITTING THE TIPS

- Screw the selected tip to the PIEZOTOME handpiece with the dynamometric wrench*;
- **Check that the spray is working properly.**

3. 2 SETTING THE MODE

a) PIEZOTOME - IMPLANT CENTER

1 being the most powerful mode

Tips	Recommended mode	Irrigation (ml/min)
LC1	1 2	60-80
LC1 90°	1 2	60-80
LC2	1 2	60-80
LC2L	1 2	60-80
LC2R	1 2	60-80
Ninja	1 2	60-80

Modes are set by the footswitch or on the front panel of the device.

b) PIEZOTOME 2 - IMPLANT CENTER 2

D1 being the most powerful mode

Tips	Recommended mode	Fine adjustment	Irrigation (ml/min)
LC1		01 02 03 04 ■■■■■ 05	60-80
LC1 90°		01 02 03 04 ■■■■■ 05	60-80
LC2		01 02 03 04 ■■■■■ 05	60-80
LC2L		01 02 03 04 ■■■■■ 05	60-80
LC2R		01 02 03 04 ■■■■■ 05	60-80
NINJA		01 02 03 04 ■■■■■ 05	60-80

Modes are set by the footswitch or on the front panel of the device.

IV - STERILIZATION

a) Lifetime

Tip shape and weight are the determining characteristics for obtaining maximal efficiency of the ultrasonic generator.

⚠ Careful monitoring of these two characteristics by the user will ensure that optimal performance of the device is maintained. Consequently, you are strongly advised not to alter the structure of the tips by filing or twisting them or modifying them in any other way.

⚠ Similarly, through normal wear, tip aging leads to alteration of its characteristics.

Always replace a tip that has been damaged by wear or by accidental impact (dropped, bent, etc.).

Replace routinely-used tips at least once to three times a year.

b) Pre-disinfection / chemical cleaning

To be carried out immediately after treatment, wearing thick gloves.

Immerse in an ultrasonic bath (enzymatic, quaternary ammonium base detergent / disinfectant solution), carefully respecting the concentration and contact time recommended by the manufacturer of the solution.

Use a solution that bears the EC marking or that is in compliance with any standard required by the national regulations.

Rinse under running water for at least 30 seconds.

c) Mechanical / chemical cleaning

To be carried out immediately after pre-disinfection / chemical cleaning, wearing thick gloves.

Scrub the products in a new bath of enzymatic, quaternary ammonium base detergent / disinfectant solution, for at least 30 seconds, and in any case until all remaining traces of contamination have been removed.

Rinse again under running water for at least 30 seconds.

d) Drying

Dry with a clean, non-woven single-use cloth, so as to remove any traces of liquid.

Pack in single-use sterilization pouches or bags in compliance with the specifications defined in

standard EN ISO 11607-1 or in any equivalent standard required by the national regulations.

e) Sterilization

The tips must be sterilized individually in an autoclave, according to the following parameters:

- autoclave: type B, compliant with standard EN 13060;
- sterilization temperature: 134°C;
- sterilization stabilization time: 18 minutes;
- pressure: 2 Bars min.

f) Storage

Then, store the sterilized products in a dry place protected from dust.

Before re-use, if the integrity of the packaging is

not certain, repackage then re-sterilize in accordance with the defined protocol.

If contamination is visible inside the bag, destroy the product.

g) Disposal of the product

Dispose of the product in receptacles for waste materials of healthcare activities involving infectious risks.

CAUTION: Not all autoclaves can reach 134°C. Not all autoclaves draw a pre-vacuum. Please refer to your autoclave manufacturer for specific sterilization instructions.

V - REPLACING TIPS

These tips must be replaced when the cutting part becomes blunt or at least once a year.

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