

Guideline: Extraction

PIEZOTOME					PIEZOTOME ²				
Maximum recommended modes					Fine setting				
	1	2	3	4	1 to 3	1 to 3	1 to 3	1 to 3	1 to 3
LC1	1	2	3	4	1 to 3	1 to 3	1 to 3	1 to 3	1 to 3
LC1 90°	1	2	3	4	1 to 3	1 to 3	1 to 3	1 to 3	1 to 3
LC2	1	2	3	4	1 to 3	1 to 3	1 to 3	1 to 3	1 to 3
LC2L	1	2	3	4	1 to 3	1 to 3	1 to 3	1 to 3	1 to 3
LC2R	1	2	3	4	1 to 3	1 to 3	1 to 3	1 to 3	1 to 3
NINJA™	1	2	3	4	1 to 3	1 to 3	1 to 3	1 to 3	1 to 3

